

Observations

| Observation | Possible Causes | Recommended Solutions |
|----------------------------|--------------------------------|---|
| Daily QC does not complete | Wrong QC bead sample | Ensure you are running SpectroFlo QC beads. |
| | Bead sample not properly mixed | Mix the bead sample. |
| | Bead sample too dilute | Concentrate the bead sample or prepare a fresh bead sample. |
| | Air bubble is sample line | Run a SIT Flush. |
| Daily QC failed | Dirty flow cell | Run a Clean Flow Cell. |



Hellmanex (let sit 15min) and MiliQ water
-then run MiliQ water for 30min on "High"
**if it still does not pass, CONTACT SRL
STAFF and wait until next business day


| Observation | Possible Causes | Recommended Solutions |
|---|-------------------------------|--|
| Data in scatter parameters appear distorted | Air bubble in flow cell | Run a SIT Flush. |
| | Air in sheath filter | Run a Purge Filter. |
| | Dirty flow cell | Run a Clean Flow Cell. |
| | Poor sample health | Check the viability of the cells. |
| High CVs | Air bubble in fluidics | Run a SIT Flush and a Purge Filter. |
| | Sample flow rate set to High | Set the sample flow rate to Low or Medium. |
| | Dirty flow cell | Run a Clean Flow Cell. If the problem persists, run a Clean Flow Cell using 25%-50% Contrad 70, followed by DI water. |
| SIT hitting bottom of well/tube | SIT Lift Distance set too low | Set the SIT Lift Distance to at least 1.5. See "Calibrating the SIT" on page 103. |

Unstable
signal/drift in
scatter signal


Dirty flow cell,
debris in fluid line,
bubbles in fluid line

Clean Flow Cell: Hellmanex (let sit 15min) and
MiliQ water
-then run MiliQ water for 30min on "High"
**if it still does not pass, CONTACT SRL STAFF
and wait until next business day

To check stability of signal: set x-axis to Time and the y-axis to which channel you had a drift in signal. The signal should remain steady with no dips or peaks



Stable



Unstable

Observations

| Observation | Possible Causes | Recommended Solutions |
|--|--|---|
| No events displayed (flow rate lower than expected) | No sample in tube | Add sample or install a new sample tube. |
| | Sample not properly mixed | Mix the sample to suspend cells/particles. |
| | Clogged SIT | Run a SIT Flush. Then run a Clean Flow Cell with 10% bleach, followed by a Clean Flow Cell with DI water. If the clog persists, replace the sample line. ★ |
| | For loaders, the SIT Lift Distance set too low (touching bottom of tube) | Increase the SIT Lift Distance. See "Calibrating the SIT" on page 103. |
| No events displayed (flow rate normal) | Insufficient gain for threshold parameter | Increase the gain for the threshold parameter. |
| | Threshold too high | Lower the threshold. |
| | Laser delay not correct | Ensure the laser delay values match those from the latest Daily QC run. See "Instrument Control" on page 44 for the laser delay location. If the values do not match, rerun Daily QC. |
| | Threshold set to incorrect parameter | Set the threshold to the appropriate parameter for the application (usually FSC). |
| | Gated plot with no data in gate | Delete or move the gate. |
| Low sample event rate | Threshold too high | Lower the threshold. |
| | Insufficient gain for threshold | Increase the gain for the threshold parameter. |
| | Sample not properly mixed | Mix the sample to suspend cells/particles. |
| | Sample too dilute | Concentrate the sample. Set the flow rate to Medium or High. |
| | Clogged SIT | Run a SIT Flush. Then run a Clean Flow Cell with 10% bleach, followed by a Clean Flow Cell with DI water. If the clog persists, replace the sample line. ★ |
| Erratic event rate | Partially blocked SIT | Run a SIT Flush. Then run a Clean Flow Cell with 10% bleach, followed by a Clean Flow Cell with DI water. |
| | Clumpy sample | Vortex, filter, or disaggregate the sample. |



★ CONTACT SRL STAFF


Concentration could also be too high. Dilute, filter and try again



This is also the cause of a DAQ board error

Loading a 96-Well Plate

- 1 If necessary, press Eject from either the QC & Setup or Acquisition module to eject the stage, bringing it forward.
- 2 Load a plate on the plate stage so that position A1 is located in the front-left corner. Insert the edge of the plate behind the metal tabs at the front of the plate stage, then press down on the back edge of the plate to secure it in the holder.



Insert back edge of plate in front of metal tabs, then press down on front edge of plate to secure.

- 3 Press Load from the QC & Setup or Acquisition module to load the plate, followed by Start to begin acquisition.

Or press Start to load the plate and begin acquisition.
- 4 To switch back to tube mode, set the Sample Delivery Mode in the Acquisition Preferences back to Tube and slide the Loader lever back.

Minimum Well Volume

When acquiring from a well, a boost is applied to deliver the sample to the flow cell as quickly as possible. This process, as well as mixing, results in the loss (consumption) of up to 75 μL of sample, as the sample is neither acquired nor recorded. We recommend a minimum volume of 150 μL of sample per well, keeping in mind that approximately 75 μL will be lost.

■ **NOTE:** With the Loader v1.5, sample consumption is reduced to <35 μL . This can be reduced further by selecting to skip the fluidics boost. For information, see “Acquisition Preferences” on page 81.

You can disable the boost to save sample. For information see “Acquisition Preferences” on page 81.