

Start-up operations (check schedules first and modify startup accordingly):

1. Start Agilent (Quanteon 1,2 and Penteon)

Daily –

- a. Check fluidics cart (NovoFlow, Waste, NovoRinse, NovoClean) and fill or empty as needed
- b. Turn on power to the cytometer (green button on left side of instrument)
- c. In NovoExpress software, go to **Instrument** tab and click **Fluidic Maintenance Sequences** → select **Startup/QC** (20min)
- d. While machines are cleaning make NovoCyte QC particles: 2 drops in 800uL PBS and 400uL of NovoRinse (working vial found in upper left pull out drawer in black fridge)
- e. After cleaning sequence is completed, click **QC test** → select the appropriate bead lot from the drop down menu
- f. After QC click **“Cleaning”**
- g. Close software and restart PC

Monthly –

- a. Instead of **Startup/QC** select **Long Clean**
- b. After QC open **“Rainbow test”** experiment. Run single peak beads to test fluidic stability
- c. In the **Instrument** tab select **“Cleaning”**
- d. Close software and restart PC

2. Start the Aurora

Daily –

- a. Fill sheath tank with MilliQ-water
- b. Empty waste container if full, add 400 mL of bleach before reattaching
- c. Replace water tube on SIT with Hellmanex 1:3
- d. Login into the SpectroFlo software and choose Acquisition from the Get Started menu
- e. Go to **Cytometer** tab and select **“Clean Flow Cell”**
 - i. Tube 1 - Hellmanex 1:3
 - ii. Wait 15 minutes before moving on to Tube #2
 - iii. Tube 2 – MilliQ-water
- f. Select **“Default”** experiment
- g. Click **“Run”** - run MilliQ-water on high for 20 minutes

- h. While water is running make tube of SpectroFlo beads: 1 drop of beads in 150uL of MilliQ-water (working vial found in red rack in butter drawer in door of black fridge)
- i. After water – go to **QC & Setup** tab
- j. Load tube of QC beads and hit **“Start”**
- k. After QC, perform **“Fluidic Shutdown”** → found under **Cytometer** tab

Monthly –

- a. Perform **“Long Clean”** under **Cytometer** tab
 - a. Replace filter with bypass tubing (found in third drawer under Aurora)
 - b. Remove tube from sheath container and take filter off the end of tube.
 - c. Place sheath tube into 10% bleach solution (2L)
 - d. Hit start in **“Long-Clean”** and follow software wizard instructions
 - e. Go through entire **“Long Clean”** with bleach solution → do not switch to water
 - f. After entire **“Long Clean”** is done. Put filter back on to end of sheath line and put back into sheath container.
 - g. Repeat **“Long Clean”** but with all steps/fluids consisting of MilliQ-water
 - h. After second **“Long Clean”** is finished, remove bypass tubing and reinstall the filter
- b. Run a tube of MilliQ water (**Default** experiment) for 20 minutes on high
- c. Perform QC per **“Daily”** protocol

3. Start the Attune

Daily –

- a. Check fluids (Waste, Focusing fluid, Wash solution, Shutdown solution)
- b. In AttuneNxT software, go to **Instrument** tab and click **Instrument Startup**
- c. After startup, load a clean dry 96-well plate into autosampler and a tube of 3mL Attune Wash solution (squeeze bottle on top of Attune) onto the SIT.
- d. Go to **Instrument** tab and select **“Quick Wash”**
- e. While machine is cleaning make Attune Performance QC particles: 3 drops in 2mL of PBS (working vial found in upper left pull out drawer in black fridge)
- f. After cleaning sequence is completed, click **Performance Test** → select the appropriate bead lot from the drop down menu
- g. After performance test load a tube of 3mL Attune Wash solution, a clean dry 96-well plate into the autosampler and click **“Quick Wash”**
- h. Close software and restart PC

4. Tyto (microfluidic sorter):

Daily –

- a. Wipe down counters with 70% EtOH
- b. Restart computer

5. Bigfoot (droplet sorter in BSL2+ suite):

If a sort is scheduled –

- a. Check bulk fluids (NERL sheath, Waste, MilliQ-water, cleaner, decontamination)
- b. Sonicate the nozzle needed for the first sort – at least 20min
- c. While nozzle is sonicating:
 - i. clean salt out of flow chamber, off charge plates, pinhole and waste catcher with MilliQ-water
 - ii. clean mirrors and back plate with Windex
 - iii. wipe down all surfaces with 70% EtOH
- d. Install sonicated nozzle and click “**Startup**” in software
- e. Let sheath run for 30 minutes then complete startup steps (droplet maintenance, stream setup, QC and drop delay)

If no sort is scheduled – every Monday, Wednesday and Friday

- f. Check bulk fluids (NERL sheath, Waste, MilliQ-water, cleaner, decontamination)
- g. Before clicking startup:
 - i. clean salt out of flow chamber, off charge plates, pinhole and waste catcher with MilliQ-water
 - ii. clean mirrors and back plate with Windex
 - iii. wipe down all surfaces with 70% EtOH
- h. click “**Startup**” in software
- i. Let sheath run for 30 minutes then complete startup steps (droplet maintenance, stream setup, QC and drop delay)
- j. Shutdown (“**Quick**”)